



**CALPROLAB**<sup>TM</sup>  
*the Calprotectin company*

## **CALPRO LAUNCHES A NEW AND IMPROVED ELISA FOR CALPROTECTIN!**

We are happy to announce that Calpro AS will shortly launch a new generation ELISA for measuring Calprotectin, our **CalproLab™** brand. Free kits are available to our customers as Research Use Only (RUO) for testing and validation. A CE-marked kit will be available for sale July 2011.

### Key features of the new ELISA:

- **25 – 2500 mg/kg** measuring range using only one dilution of faecal extracts
- Simplified extraction procedure using Faecal extraction device
- More user friendly ELISA: shorter assay time and the same volume for all steps in the assay
- Six standards and two positive controls
- Monoclonal – polyclonal antibody combination
- Validated by leading European clinics and laboratories
- The wide assay range makes the CalproLab™ kit ideal for both treatment regimen monitoring, relapse prediction and discrimination between IBS and IBD
- In addition to faecal extracts, the assay is also validated for testing of **plasma/serum samples**

### Additional notes:

- Validation of our kits is important – we welcome requests for kits to participate in clinical studies!
- The new ELISA will have product number CALP0170 (ALP) and CALP0270 (HRP) – note that our current kits, CAL0100 (ALP) and CAL0300 (HRP), will still be available

Please contact Janne K Eidem [jeidem@calpro.no](mailto:jeidem@calpro.no) for technical information and Magdalena Urbanik [murbanik@calpro.no](mailto:murbanik@calpro.no) for commercial questions.

Steinar Ask Magnussen, CEO (sign.)

Lysaker 23.05.11

CALPRO AS

